

PROCEDURE TO FOLLOW FOR THE AMNIOTIC FLUID COLLECTION

SAMPLE

- Collect 1 ml per week of amenorrhoea on sterile vials (available on request)
- Store at room temperature (<40°C)
- Identify the vial(s)

ADMINISTRATIVE

- Please include with the sample:
 - 1) The duly completed medical prescription form
 - 2) The signed certificate from the prescribing physician and the signed consent form from the patient
 - 3) Clinical pathology results report or sonographic report according to the clinical indications.